

# EVALUATION OF STEMULATE™ FOR THE EXPANSION AND CRYOPRESERVATION OF HUMAN MESENCHYMAL STEM CELLS

Thomas Heathman<sup>1</sup>, Karen Coopman<sup>1</sup>, Alvin Nienow<sup>1,2</sup>, Bo Kara<sup>3</sup>, Christopher Hewitt<sup>1,4</sup>


<sup>1</sup>Centre for Biological Engineering, Loughborough University, UK  
<sup>3</sup>FUJIFILM Diosynth Biotechnologies, Billingham, UK

<sup>2</sup>School of Chemical Engineering, University of Birmingham, UK  
<sup>4</sup>School of Life and Health Sciences, Aston University, UK

## Overview


The cost-effective production of AD- & BM-hMSC therapies will require scalable processes with the ability to consistently deliver hMSCs to the clinic. In addition to lot-to-lot variability, there are additional process constraints on the use of bovine-derived serum such as limited supply, potential for pathogen transmission and increasing cost at scale. These considerations are likely to limit the scalability of hMSC manufacturing processes based on bovine-serum and therefore must be addressed. Stemulate™ (human platelet lysate) has been proposed as a viable alternative and this work aims to evaluate this potential.

## MONOLAYER




- Adipose-derived (AD-) and bone marrow-derived (BM-) hMSCs have been successfully expanded in DMEM supplemented with 10% bovine serum (FBS) and 5% human platelet lysate (HPL) over six passages (Figure A & B).
- Despite using half the amount of HPL, AD-hMSC expansion was significantly increased compared with FBS (Figure A & B).
- Expression of VEGF has been implicated in the promotion of angiogenesis and has remained present (Figure C).
- CCL2 is an important marker for hMSC immune modulation and has been retained throughout expansion (Figure D).
- Per cell metabolic flux of BM-hMSCs suggests utilisation of similar pathways in FBS and HPL-based culture (Figure E & F).

## MICROCARRIER




- BM-hMSCs were expanded on non-porous, plastic microcarriers in 100 mL spinner flasks over six days in DMEM supplemented with FBS or HPL.
- Process transfer to a scalable platform is an important step in the development of cost-effective processes.
- HPL increased the expansion rate of both hMSC donors, eliminating the lag phase observed in FBS (Figure G).
- Images show BM-hMSC growth on microcarriers in suspension (Figure H). Live cells stained with Calcein AM fluorophore (Figure I).
- BM-hMSC retained all desired characteristics post microcarrier harvest including CFU potential (Figure J).




Process map for the expansion, harvest, downstream processing and preservation of hMSCs from monolayer & microcarrier processes.

- Post-microcarrier and monolayer characterisation of hMSCs cultured in FBS and HPL-based medium demonstrate the maintenance of key characteristics according to the ISCT minimum criteria, namely:

- Adherence to tissue culture plastic (Figure S & T).
- Tri-lineage differentiation down osteo-, adipo- & chondrogenic lineages (Figure R).
- Co-expression of positive immunophenotype markers CD73, 90, 105 & 44 (Figure Q).
- Co-expression of negative immunophenotype markers CD45, 34, 11b, 19 & HLA-DR (Figure Q).



## ISCT CRITERIA



## PRESERVATION

## Conclusion & Perspective

Demonstrating the amenability of Stemulate™ as a replacement for bovine serum in the expansion and cryopreservation of hMSCs provides a step forward in developing scalable manufacturing processes. The increased growth kinetics and batch consistency associated with Stemulate™ combined with the reduced volume required for culture make it a viable alternative to bovine serum for future process development.

## Contact

t.r.j.heathman@lboro.ac.uk  
 +44(0)1509 564891  
 @Thomas\_Heathman

