

Keep a Gauge on the Pressure



Compass[®] CT

DISPOSABLE PRESSURE TRANSDUCER

The Compass CT is a disposable pressure transducer featuring an integrated digital display. It provides physicians more control by monitoring physiological pressure and displaying pressure during the infusion of therapeutic agents. The Compass CT Pressure Transducer provides:

- **Precision:** measures pressures from -199 to +999 mm Hg
- **Accuracy:** displays pressure $\pm 3\%$ above 50 mm Hg
- **Consistency:** provides zero drift ± 1 mm Hg per 4 hours
- **Versatility:** allows for multiple readings throughout one procedure

RegenDevices

Tools for therapeutic delivery

RegenDevices are tools for delivering therapeutic agents to the target, designed for precision, accuracy, and consistency.

Customer Service

Phone: +1 317.917.3450

Toll-free in the United States: +1 800.265.0945

E-mail: customer-service@cookregentec.com

Cook Regentec

1102 Indiana Avenue
Indianapolis, IN 46202 U.S.A.

This product may not be available in all markets.

Product	Pressure Range mm Hg	RPN	Order Number
Compass [®] CT Disposable Pressure Transducer	-199 to +999	CCT001	G29279



www.cookregentec.com

Keep a Gauge on the Pressure



Compass® CT Port

DISPOSABLE PRESSURE TRANSDUCER

The Compass CT Port is a disposable pressure transducer featuring an integrated digital display. It provides physicians more control by monitoring physiological pressure and displaying pressure during the infusion of therapeutic agents. Its port allows for passing of wire guides and other devices. The Compass CT Port Pressure Transducer provides:

- **Precision:** measures pressures from -199 to +999 mm Hg
- **Accuracy:** displays pressure $\pm 3\%$ above 50 mm Hg
- **Consistency:** provides zero drift ± 1 mm Hg per 4 hours
- **Versatility:** allows for passing wire guides or devices up to 0.038 inch (0.97 mm) in diameter

Product	Pressure Range mm Hg	RPN	Order Number
Compass® CT Port Disposable Pressure Transducer	-199 to +999	CCTP001	G29278

RegenDevices

Tools for therapeutic delivery

RegenDevices are tools for delivering therapeutic agents to the target, designed for precision, accuracy, and consistency.

Customer Service

Phone: +1 317.917.3450

Toll-free in the United States: +1 800.265.0945

E-mail: customer-service@cookregentec.com

Cook Regentec

1102 Indiana Avenue
Indianapolis, IN 46202 U.S.A.

This product may not be available in all markets.



www.cookregentec.com