

Keep a Gauge on the Pressure



Compass[®] CT

DISPOSABLE PRESSURE TRANSDUCER

The Compass CT is a disposable pressure transducer featuring an integrated digital display. It provides physicians more control by monitoring physiological pressure and displaying pressure during the infusion of therapeutic agents. The Compass CT Pressure Transducer provides:

- **Precision:** measures pressures from -199 to +999 mm Hg
- **Accuracy:** displays pressure $\pm 3\%$ above 50 mm Hg
- **Consistency:** provides zero drift ± 1 mm Hg per 4 hours
- **Versatility:** allows for multiple readings throughout one procedure

Product	Pressure Range mm Hg	RPN	Order Number
Compass [®] CT Disposable Pressure Transducer	-199 to +999	CCT001	G35061

RegenDevices

Tools for therapeutic delivery

RegenDevices are tools for delivering therapeutic agents to the target, designed for precision, accuracy, and consistency.

Customer Service

United States & Canada
Toll Free: +1.800.457.4500
Phone: +1.812.339.2235
CustomerSupport@CookMedical.com

Beyond North America
Phone: +1.317.917.3450
Customer-Service@CookRegentec.com

Cook Regentec
1102 Indiana Avenue
Indianapolis, IN 46202 U.S.A.

This product may not be available in all markets.



Keep a Gauge on the Pressure



Compass® CT Port

DISPOSABLE PRESSURE TRANSDUCER

The Compass CT Port is a disposable pressure transducer featuring an integrated digital display. It provides physicians more control by monitoring physiological pressure and displaying pressure during the infusion of therapeutic agents. Its port allows for passing of wire guides and other devices. The Compass CT Port Pressure Transducer provides:

- **Precision:** measures pressures from -199 to +999 mm Hg
- **Accuracy:** displays pressure $\pm 3\%$ above 50 mm Hg
- **Consistency:** provides zero drift ± 1 mm Hg per 4 hours
- **Versatility:** allows for passing wire guides or devices up to 0.038 inch (0.97 mm) in diameter

Product	Pressure Range mm Hg	RPN	Order Number
Compass® CT Port Disposable Pressure Transducer	-199 to +999	CCTP001	G35060

RegenDevices

Tools for therapeutic delivery

RegenDevices are tools for delivering therapeutic agents to the target, designed for precision, accuracy, and consistency.

Customer Service

United States & Canada
Toll Free: +1.800.457.4500
Phone: +1.812.339.2235
CustomerSupport@CookMedical.com

Beyond North America
Phone: +1.317.917.3450
Customer-Service@CookRegentec.com

Cook Regentec
1102 Indiana Avenue
Indianapolis, IN 46202 U.S.A.

This product may not be available in all markets.



cookregentec.com